

THE MOST STABLE CUSHION FOR WHEELCHAIR USERS

Ischial well - Two compartments
for superior pressure reduction in the IT-coccyx area & extra positioning possibilities

Trochanter support
firmer density to provide lateral stability | redistributes pressure to help create an ischial well



Front compartments provide leg alignment and stability | easily adjustable for individual positioning needs

Mid height zip closure for easy adjustment when client is seated | cushion keeps its shape, enhancing the user's skin protection

IMPROVED DESIGN



Expanding the user's reach

Nine adjustable compartments for maximum stability

Extra positioning of the thighs

Sometimes it is necessary to reduce pressure in the vulnerable area around the IT's (sitting bones) and coccyx (tailbone). The Vector 02 does just that. Its unique compartment layout further provides extra pelvic stability, excellent leg guidance, optimal skin protection and comfort.

100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation. The Vector 02 offers support where needed, resulting in improved functionality of the upper body and expanding reach.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus	▲▲▲▲▲
Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲△
Micro climate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲▲
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲▲
Allows easy transfers	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲△
Low maintenance level	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Low cushion weight	▲▲▲▲▲
Available size range	▲▲▲▲▲
Maximum user weight	▲▲▲▲▲
Low seat to floor height	▲▲▲▲▲
Auto adjustments	▲▲▲▲△

USED WITH


Spinal cord injuries, neuromuscular diseases, hypermobility syndrome, Ehlers Danlos Syndrome (EDS)

TECHNICAL SPECIFICATIONS

Height	6 and 10 cm
Average weight	725 grams
Max. user weight	250 kg
Machine washable	60°C max.

CUSHION SIZES

■ Standard (Sizes >53 cm height 10 only)
■ On request



NEW



VICAIR[®]


O2

The Vicair O2 cushions tackle the hygiene issues caused by incontinence. The open structure of the O2 cushions allows the cushions to be easily washed in any washing machine. Without the need for a laundry bag, and without removing the SmartCells. These cushions are ideal for both institutional and home use, in situations where incontinence, temperature and humidity is an issue.




100% MACHINE WASHABLE

The cushion is, with the SmartCells inside, fully machine washable at 60°C



100% BREATHABLE

The open structure of the cushion provides optimal temperature and moisture regulation



SOLID SHAPE

This new cushion model keeps its shape, enhancing the user's skin protection.



WIDE RANGE

Available in two heights; 145 standard sizes and custom sizes available on request

Fully machine washable for optimal hygiene

100% breathable for perfect temperature and moisture regulation

Ready for use



relief by air