

MAATTABEL ANATOMIC BACK KUSSEN



Beschikbare maten:

		BREEDTE (BB)				
inch		12	14	16	18	20
cm		30	35	40	45	50
HOOGTE (HH)	12¾					
	16					
	19					
	22½					

 Standaard maat  Speciale maat

Uit voorraad leverbaar

4 weken levertijd

MAAT = BBHH (breedte x hoogte)

Strap Back

Art.code	Min-max afstand tussen rugbuizen van de rolstoel		Dikte Strap Back	
	cm	inches	cm	inches
SB/0132	30	12	32	12¾
SB/0232	35-40	14-16	32	12¾
SB/0332	45-50	18-20	32	12¾
SB/0140	30	12	40	16
SB/0240	35-40	14-16	40	16
SB/0340	45-50	18-20	40	16
SB/0148	30	12	48	19
SB/0248	35-40	14-16	48	19
SB/0348	45-50	18-20	48	19
SB/0156	30	12	56	22½
SB/0256	35-40	14-16	56	22½
SB/0356	45-50	18-20	56	22½

