

PERFECT TO USE IN COMBINATION WITH LIBERTY BACK AND ANATOMIC BACK

Hook and look strips to attach cushion

Adjustable straps to adjust to back contours. Number of straps depends on Strap back length

Easy and secure full length hook and loop attachment to wheelchair canes, with screw fixation for additional security

Double zip at the front for addition of lateral support

Double zip at the back to allow fixation of stabilisation bar



- ▲ **Easy to apply on manual wheelchairs**
- ▲ **Perfect positioning and comfort qualities**
- ▲ **Very lightweight**

The Vicair Strap back is a perfect base for Vicair wheelchair back cushions. Quickly installed and combined with your choice of Vicair back cushion it is a high quality alternative to your standard wheelchair back system. The adjustable straps enhance the function of the Vicair back cushions, by offering thoracic stability and backrest positioning.

USED WITH

Vicair Liberty back, Vicair Anatomic back, Vicair Multifunctional O2 back (in combination with the Vicair Back Connector)

TECHNICAL SPECIFICATIONS

Max. user weight 150 kg