

**EXCELLENT FOR MANUAL WHEELCHAIR USERS WITH BACK PAIN**



**Air technology** inside for comfort, support and pressure redistribution



Combine with easy-install adjustable Vicair Strap back for positioning and balance



- ▲ **Easy to apply on a standard (folding) wheelchair**
- ▲ **Automatically adapts to change, including reclining back support**
- ▲ **Excellent for users with back pain**

The Vicair Liberty back cushion changes your wheelchair back from standard into extremely comfortable. Combined with the easy-install Vicair Strap back, the cushion offers perfect pain relief, and follows and gently supports postural and/or anatomical anomalies. This extremely lightweight back cushion is easy in transport, as the wheelchair can still be folded when the cushion is attached to it. An easy solution to create comfort and support.

**PRODUCT PERFORMANCE**

**Skin protection / Anti decubitus**

Pressure redistribution	▲▲▲△△
Shear force reduction	▲▲▲△△
Microclimate control	▲▲△△△
Impact & vibration damping	▲▲▲△△

**Positioning**

Lateral stability	▲▲▲△△
Frontal stability	N/A
Sliding reduction	▲▲▲△△
Allows manual adjustments	N/A
Increasing endurance	▲▲▲△△

**Maintenance & Safety**

Durability	▲▲▲△△
Reliability	▲▲▲△△
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲△△
Level of hygiene	▲▲▲△△

**Functionality & Application**

Comfort	▲▲▲△△
Ultra light weight	▲▲▲△△
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲△△△
Maximum user weight	▲▲▲△△

**TECHNICAL SPECIFICATIONS**

Length	45 cm
Average weight	350 grams
Max. user weight	150 kg
Loss of seath depth	3 cm

**CUSHION SIZES**

Available in 3 sizes. Mentioned cm/inches are the min-max distance between the back canes of the wheelchair.

	cm	inches
LBS/1	36-40	14¼ - 16
LBS/2	40-44	16 - 17½
LBS/3	44-48	17½ - 19

**AVAILABLE COVERS**

Comfair Cover.