


EXTREMELY USER-FRIENDLY AND HIGH PERFORMANCE SKIN PROTECTION




Colour coded label
Easily shows if cushion and cover are the same size

Safe to use in all directions:
top = bottom; left = right;
front = back.



4 equal compartments



- ▲ High performance skin protection by even pressure re-distribution
- ▲ Extreme ease of use
- ▲ Supporting the pelvis and thighs for a stable sitting posture

The Vicair 4 consists of 4 identical individual compartments designed to contour under loading. The cushion is always placed correctly, as all four sides can be either front or back even the top and bottom are identical. The configuration of the four compartments is designed to work together to strike a unique balance between comfort, preserving skin and pelvic support. The four identical individual compartments allow for a foolproof correct positioning on the seating surface. Vicair 4 cushions and covers have the unique feature of colour coded labels, easily showing if the cushion and cover are the same size. This is considered to be beneficial in environments where many wheelchair cushions are in use at once, such as care centres.

AVAILABLE COVERS

Comfair Cover and Incotec Cover.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲△△
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	N/A
Increasing endurance	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲▲
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Ultra light weight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲△△
Maximum user weight	▲▲▲▲△

TECHNICAL SPECIFICATIONS

Thickness	8 cm
Average weight	800 grams
Max. user weight	250 kg
Machine washable	60°C max.

CUSHION SIZES

- Available
- Not available

