

THE BEST BACK SUPPORT FOR USERS OF MANUAL TILT-IN-SPACE WHEELCHAIRS AND POWERCHAIRS



- ▲ Supporting different postures
- ▲ Follows the human spine alignment
- ▲ Optimal pressure redistribution

Excellent pressure redistribution, comfort and support make this back support cushion perfect in a large variety of situations. It can be used in several types of back systems and in tilt-in-space chairs with adjustable seat depth. The cushion is easily fixed onto the solid back support by use of hook and loop strips. With the Vicair Back Connector the Multifunctional O2 Back can also be used on tension adjustable backs such as the Vicair Strap Back. The cushion is developed to support an upright position by creating extra support at lumbar thoracic extension of the spine. It also automatically adapts to the body's contours. The cushion can be adjusted to accommodate optimal positioning and/or pressure management, this while the user is still leaning against the cushion.

AVAILABLE COVERS
Comfair Cover.

| PRODUCT PERFORMANCE | |
|---|-------|
| Skin protection / Anti decubitus | |
| Pressure redistribution | ▲▲▲▲▲ |
| Shear force reduction | ▲▲▲▲△ |
| Microclimate control | ▲▲▲▲▲ |
| Impact & vibration damping | ▲▲▲▲▲ |
| Positioning | |
| Lateral stability | ▲▲▲▲△ |
| Frontal stability | N/A |
| Sliding reduction | ▲▲▲▲▲ |
| Allows manual adjustments | ▲▲▲▲▲ |
| Increasing endurance | ▲▲▲▲▲ |
| Maintenance & Safety | |
| Durability | ▲▲▲▲▲ |
| Reliability | ▲▲▲▲▲ |
| Functionality when misused | ▲▲▲▲△ |
| Easy to maintain | ▲▲▲▲△ |
| Level of hygiene | ▲▲▲▲▲ |
| Functionality & Application | |
| Comfort | ▲▲▲▲▲ |
| Ultra lightweight | ▲▲▲▲△ |
| Allows easy transfers | ▲▲▲▲△ |
| Available standard sizes | ▲▲▲▲△ |
| Maximum user weight | ▲▲▲▲△ |

| TECHNICAL SPECIFICATIONS | |
|--------------------------|-----------|
| Length | 36-60 cm |
| Average weight | 975 grams |
| Max. user weight | 150 kg |
| Machine washable | 60°C max. |
| Loss of seat depth | 3 cm |

CUSHION SIZES
■ Standard
■ On request

