

LIGHTWEIGHT POSITIONING SUPPORT, COMFORTABLE AND EASY TO USE



- ▲ More stability thanks to the shape of the integrated foam base
- ▲ Air valves automatically adapt the cushion to your body shape and weight
- ▲ Comfortable and rain proof

The Vicair Liberty Profile is an easy to use and comfortable wheelchair cushion. Designed for comfort in combination with a low to medium level of support and pressure injury prevention. The Vicair Liberty Profile has an integrated foam base with a semi-circle cut out at the location of the sitting bones. The combination of SmartCells and foam-base supports pelvic stability and offers a high level of comfort in the IT-Coccyx zone. The added second valve allows an even easier automatic adaption to the body shape. This Vicair cushion is ready-for-use, rain proof and very lightweight.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲△△
Shear force reduction	▲▲△△△
Microclimate control	▲▲△△△
Impact & vibration damping	▲▲▲△△

Positioning

Lateral stability	▲▲▲△△
Frontal stability	▲▲▲△△
Sliding reduction	▲▲▲△△
Allows manual adjustments	N/A
Increasing endurance	▲▲▲△△

Maintenance & Safety

Durability	▲▲▲△△
Reliability	▲▲▲△△
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲△△
Level of hygiene	▲▲▲△△

Functionality & Application

Comfort	▲▲▲△△
Ultra lightweight	▲▲▲△△
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲△△
Maximum user weight	▲▲▲△△

TECHNICAL SPECIFICATIONS

Thickness	6 cm
Average weight	600 grams
Max. user weight	150 kg

CUSHION SIZES

- Standard
- Not available

Inch	width										
	14	15	16	17	18	19	20	21	22	23	24
CM	35	38	40	43	45	48	50	53	55	58	60
14	35										
15	38										
16	40										
17	43										
18	45										
19	48										
20	50										
21	53										
22	55										
23	58										
24	60										

AVAILABLE COVERS

Air Exchange Cover and Incotec Cover.