



PRESSURE RELIEF AND STABILITY WITH A DETACHABLE MEDIAL THIGH SUPPORT

Follows the body's contours for equal pressure distribution | support IT-coccyx area | managing pelvic asymmetry



Medial thigh support

Thigh support



- ▲ Ideal for users with adducted conditions
- ▲ Detachable medial thigh support
- ▲ Cushion and medial thigh support covered separately for ease of use

The Vicair Pommel cushion O2 consists of 5 individual compartments including a detachable medial thigh support. This cushion is ideal for individuals with adducted conditions, where the knees or thighs are pulled together. The detachable medial thigh support gently separates the knees and/or thighs relieving the area of pressure that could normally result in pain. The SmartCell filling ensures a softer padded surface between the knees and/or thighs, redistributing pressure in this area and therefore also reducing the risk for pressure injuries. The detachable medial thigh support can be removed to facilitate forward transfers or catheterisation whilst remaining in the wheelchair.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲△
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲△
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲△
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲▲△
Level of hygiene	▲▲▲▲▲

Functionality & Application


Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲▲▲
Maximum user weight	▲▲▲▲△

TECHNICAL SPECIFICATIONS

Thickness	10 cm
Average weight	850 grams
Max. user weight	250 kg
Machine washable	60°C max.

CUSHION SIZES

- Standard
- On request



AVAILABLE COVERS

Comfair Cover and Incotec Cover.