


**EASY TO USE AND COMFORTABLE**

Front and rear compartments | reduces tendency to slide

High quality pressure redistribution | Immersion



- ▲ Two compartment layout limits sliding
- ▲ Comfortable
- ▲ Optimal immersion and envelopment

The Vicair Twin O2 is an easy to use and comfortable wheelchair cushion. It is designed for those who require comfort in combination with a high level of support and pressure injury prevention.

The two compartments of the cushion provide an even pressure distribution and help to prevent sliding forward.

100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation.

The Vicair Twin O2 is lightweight and ready-for-use.

**PRODUCT PERFORMANCE**

**Skin protection / Anti decubitus**

Pressure redistribution	▲▲▲▲△
Shear force reduction	▲▲▲△△
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲△

**Positioning**

Lateral stability	▲▲▲△△
Frontal stability	▲▲▲▲△
Sliding reduction	▲▲▲△△
Allows manual adjustments	▲▲▲△△
Increasing endurance	▲▲▲△△

**Maintenance & Safety**

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

**Functionality & Application**

Comfort	▲▲▲▲△
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲△△△
Maximum user weight	▲▲▲▲▲

**TECHNICAL SPECIFICATIONS**

Thickness	6 and 10 cm
Average weight	650 - 935 grams
Max. user weight	250 kg
Machine washable	60°C max.

**CUSHION SIZES**

- Standard
- Special Size

Inch	width										
	14	15	16	17	18	19	20	21	22	23	24
CM	35	38	40	43	45	48	50	53	55	58	60
14	35										
15	38										
16	40										
17	43										
18	45										
19	48										
20	50										
21	53										
22	55										
23	58										
24	60										

Displayed sizes are for the 10cm high model. For 6cm high model contact dealer.

070069-OTW-EN-REV01/2024

**AVAILABLE COVERS**

Vicair O2 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.