

ESPECIALLY DEVELOPED FOR BARIATRICS

Follows the body's contours
for equal pressure distribution | support IT-coccyx area | managing pelvic asymmetry

Pre-ischial bar limits
sliding forward | converts shear into normal forces



Thigh support

PRESSURE INJURY TREATMENT PLAN

100% BREATHABLE

SOLID SHAPE

MACHINE WASHABLE

READY FOR USE

- ▲ **User weight up to 500kg**
- ▲ **Shaped like the Vicair Adjuster 02**
- ▲ **Wider size range than other Vicair products**

The Vicair XXtra 02 consists of 5 individual compartments, designed to contour under loading and contains especially designed – large – SmartCells to be able to support bariatric individuals. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials, coccyx and trochanters, positioning the pelvis and supporting the femurs. The wheelchair cushion is able to support a body weight up to 500KG. For specific needs the filling of each of the 5 compartments can be easily adjusted, allowing you to change the density of each compartment.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲▲
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲▲△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲▲△
Maximum user weight	▲▲▲▲▲

TECHNICAL SPECIFICATIONS

Thickness	10 cm
Average weight	900 grams
Max. user weight	500 kg
Machine washable	60°C max.

CUSHION SIZES

- Standard
- On request

Inch	width					
	22	24	26	28	30	
CM	55	60	65	70	75	
18	45					
20	50					
22	55					
24	60					
26	65					
28	70					
30	75					

AVAILABLE COVERS

Vicair 02 Top Cover, Comfair Cover and Incotec Cover.