

OFFLOADING PRESSURE IN THE HIGH RISK AREA

PRESSURE INJURY TREATMENT PLAN
100% BREATHABLE
SOLID SHAPE
MACHINE WASHABLE
READY FOR USE

- ▲ For completely offloading of the rectal area
- ▲ Adjustable non-contact area
- ▲ Stable positioning

The Vicair Centre Relief 02 consists of 5 individual compartments, designed to contour under loading. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials and trochanters, positioning the pelvis and supporting the femurs. The unique feature is the circular compartment towards the rear of the cushion. Designed especially for (partial) offloading of the coccyx/rectal area. For specific needs the filling of each of the 5 compartments can be easily adjusted, for individuals who: (1) require (partial) offloading of the rectal area, (2) have coccyx problems, (3) have localized deep pressure injuries.

AVAILABLE COVERS

Vicair 02 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.

PRODUCT PERFORMANCE
Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲▲
Shear force reduction	▲▲▲▲△
Microclimate control	▲▲▲▲▲
Impact & vibration damping	▲▲▲▲▲

Positioning

Lateral stability	▲▲▲▲△
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲△△
Allows manual adjustments	▲▲▲▲▲
Increasing endurance	▲▲▲▲△

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲▲▲
Level of hygiene	▲▲▲▲▲

Functionality & Application

Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲▲
Allows easy transfers	▲▲▲△△
Available standard sizes	▲▲▲△△
Maximum user weight	▲▲▲▲△

TECHNICAL SPECIFICATIONS

Thickness	10 cm
Average weight	750 grams
Max. user weight	250 kg
Machine washable	60°C max.

CUSHION SIZES

- Standard
- On request

		width												
Inch	CM	14	15	16	17	18	19	20	21	22	23	24		
		35	38	40	43	45	48	50	53	55	58	60		
14	35													
15	38													
16	40													
17	43													
18	45													
19	48													
20	50													
21	53													
22	55													
23	58													
24	60													

070069-OCR-EN-REV01/2025